

— ABLYBIO, Help Your Research



PELP1 (YD32832) Rabbit mAb

货号: **AYD11385**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	120kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,MCF7,SKOV3,NIH/3T3,Rat testis
细胞定位	Nucleus, nucleolus, nucleoplasm, Cytoplasm
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a transcription factor which coactivates transcription of estrogen receptor responsive genes and corepresses genes activated by other hormone receptors or sequence-specific transcription factors. Expression of this gene is regulated by both members of the estrogen receptor family. This gene may be involved in the progression of several types of cancer. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013]
基因ID	27043
基因名	PELP1
Swiss	Q8IZL8 (https://www.uniprot.org/uniprotkb/Q8IZL8/entry)
别名	PELP1 (YD32832),PELP1 (YD32832) Rabbit mAb,PELP1,Modulator of non-genomic activity of estrogen receptor,Transcription factor HMX3,HMX3,MNAR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)