

— ABLYBIO, Help Your Research



# MSH2 (YD12019) Rabbit mAb

货号: **AYD11396**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	105kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	HeLa, A-549
细胞定位	Nucleus, Chromosome
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	This locus is frequently mutated in hereditary nonpolyposis colon cancer (HNPCC). When cloned, it was discovered to be a human homolog of the E. coli mismatch repair gene mutS, consistent with the characteristic alterations in microsatellite sequences (RER+ phenotype) found in HNPCC. Two transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4436
基因名	MSH2
Swiss	P43246 ( <a href="https://www.uniprot.org/uniprotkb/P43246/entry">https://www.uniprot.org/uniprotkb/P43246/entry</a> )
别名	MSH2 (YD12019),MSH2 (YD12019) Rabbit mAb,MSH2,MutS protein homolog 2

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)