

— ABLYBIO, Help Your Research



CDC37 (YD31139) Rabbit mAb

货号: **AYD11420**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IP |
| 推荐浓度 | |
| 理论分子量 | 44kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | Raji,HeLa,HepG2,BT474,HT-1080 |
| 细胞定位 | Cytoplasm |
| 纯化 | 亲和纯化 |



抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | The protein encoded by this gene is highly similar to Cdc 37, a cell division cycle control protein of <i>Saccharomyces cerevisiae</i> . This protein is a molecular chaperone with specific function in cell signal transduction. It has been shown to form complex with Hsp90 and a variety of protein kinases including CDK4, CDK6, SRC, RAF-1, MOK, as well as eIF2 alpha kinases. It is thought to play a critical role in directing Hsp90 to its target kinases. |
| 基因ID | 11140 |
| 基因名 | CDC37 |
| Swiss | Q16543 (https://www.uniprot.org/uniprotkb/Q16543/entry) |
| 别名 | CDC37 (YD31139), CDC37 (YD31139) Rabbit mAb, CDC37, Hsp90 chaperone protein kinase-targeting subunit, p50Cdc37, CDC37A |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)