

— ABLYBIO, Help Your Research



# Surfactant protein D (YD15785) Rabbit mAb

货号: **AYD11433**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	38kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted, extracellular space, extracellular matrix, surface film
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Contributes to the lung's defense against inhaled microorganisms, organic antigens and toxins. Interacts with compounds such as bacterial lipopolysaccharides, oligosaccharides and fatty acids and modulates leukocyte action in immune response. May participate in the extracellular reorganization or turnover of pulmonary surfactant. Binds strongly maltose residues and to a lesser extent other alpha-glucosyl moieties
------	--

基因ID	3134
基因名	Sftpd
Swiss	P50404 ( <a href="https://www.uniprot.org/uniprotkb/P50404/entry">https://www.uniprot.org/uniprotkb/P50404/entry</a> )
别名	Surfactant protein D (YD15785), Surfactant protein D (YD15785) Rabbit mAb, Sftpd, Lung surfactant protein D, Sftp4

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)