

— ABLYBIO, Help Your Research



# RAP80 (YD33241) Rabbit mAb

货号: **AYD11498**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	80kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,293T,Mouse thymus,Rat brain
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a nuclear protein that interacts with Brca1 (breast cancer 1) in a complex to recognize and repair DNA lesions. This protein binds ubiquitinated lysine 63 of histone H2A and H2AX. This protein may also function as a repressor of transcription. Alternative splicing results in multiple transcript variant s.
------	--

基因ID	51720
基因名	UIMC1
Swiss	Q96RL1 ( <a href="https://www.uniprot.org/uniprotkb/Q96RL1/entry">https://www.uniprot.org/uniprotkb/Q96RL1/entry</a> )
别名	RAP80 (YD33241),RAP80 (YD33241) Rabbit mAb,UIMC1,Receptor-associated protein 80,Retinoid X receptor-interacting protein 110,Ubiquitin interaction motif-containing protein 1,RAP80,RXRIP110

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)