

— ABLYBIO, Help Your Research



PCNA (YD14426) Rabbit mAb

货号: **AYD11523**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	29kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,MCF7,293T,COS-7,Rat testis,C2C12
细胞定位	Nucleus
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is found in the nucleus and is a cofactor of DNA polymerase delta. The encoded protein acts as a homotrimer and helps increase the processivity of leading strand synthesis during DNA replication. In response to DNA damage, this protein is ubiquitinated and is involved in the RAD6-dependent DNA repair pathway. Two transcript variants encoding the same protein have been found for this gene. Pseudogenes of this gene have been described on chromosome 4 and on the X chromosome.
基因ID	5111
基因名	PCNA
Swiss	P12004 (https://www.uniprot.org/uniprotkb/P12004/entry)
别名	PCNA (YD14426),PCNA (YD14426) Rabbit mAb,PCNA,Cyclin,Proliferating cell nuclear antigen

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)