

— ABLYBIO, Help Your Research



# WAVE 1 (YD16838) Rabbit mAb

货号: **AYD11567**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	62kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,DU145,HeLa,K-562,293T
细胞定位	Cytoplasm, cytoskeleton, Synapse, Cell junction, focal adhesion
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene, a member of the Wiskott-Aldrich syndrome protein (WASP)-family, plays a critical role downstream of Rac, a Rho-family small GTPase, in regulating the actin cytoskeleton required for membrane ruffling. It has been shown to associate with an actin nucleation core Arp2/3 complex while enhancing actin polymerization in vitro. Wiskott-Aldrich syndrome is a disease of the immune system, likely due to defects in regulation of actin cytoskeleton. Multiple alternatively spliced transcript variants encoding the same protein have been found for this gene.
基因ID	8936
基因名	WASF1
Swiss	Q92558 ( <a href="https://www.uniprot.org/uniprotkb/Q92558/entry">https://www.uniprot.org/uniprotkb/Q92558/entry</a> )
别名	WAVE 1 (YD16838), WAVE 1 (YD16838) Rabbit mAb, WASF1, Protein WAVE-1, Verprolin homology domain-containing protein 1, Wiskott-Aldrich syndrome protein family member 1, KIAA0269, SCAR1, WAVE1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)