

VPS28 (YD34410) Rabbit mAb

货号: **AYD11578**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	25kDa/25kDa/26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane, Late endosome membrane, Cytoplasm, Endosome
纯化	

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Component of the ESCRT-I complex, a regulator of vesicular trafficking process Component of the ESCRT-I complex, a regulator of vesicular trafficking process Component of the endosomal sorting complex required for transport-I (ESCRT-I) complex, a regulator of vesicular trafficking process. The ESCRT-I complex participates in the sorting of membrane proteins into vesicles that bud into late endosomal compartments called multivesicular bodies, and is also involved in the budding of many viruses
基因ID	51160
基因名	VPS28, Vps28
Swiss	Q9UK41, Q9D1C8, B5DEN9

别名	VPS28 (YD34410)
----	-----------------

产品验证

Jurkat
kDa
250 -
150 -
100 -
75 -
50 -
37 -
25 - ←
20 -
15 -
10 -

实验步骤

访问官网浏览详情: www.ablybio.cn