

— ABLYBIO, Help Your Research



Matrin 3 (YD11375) Rabbit mAb

货号: **AYD11647**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	95kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,MCF7,K-562,NIH/3T3
细胞定位	Nucleus matrix
纯化	亲和纯化




抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a nuclear matrix protein, which is proposed to stabilize certain messenger RNA species. Mutations of this gene are associated with distal myopathy 2, which often includes vocal cord and pharyngeal weakness. Alternatively spliced transcript variants, including read-through transcripts composed of the upstream small nucleolar RNA host gene 4 (non-protein coding) and matrin 3 gene sequence, have been identified. Pseudogenes of this gene are located on chromosomes 1 and X.
基因ID	9782
基因名	MATR3
Swiss	P43243 (https://www.uniprot.org/uniprotkb/P43243/entry)
别名	Matrin 3 (YD11375),Matrin 3 (YD11375) Rabbit mAb,MATR3,KIAA0723

产品验证





实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)