

FGL1 (YD35877) Rabbit mAb

货号: **AYD11650**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC IP
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2
细胞定位	Secreted
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	Fibrinogen-like protein 1 (FGL1), a liver-secreted protein, is a major LAG-3 functional ligand independent from MHC-II. The FGL1-LAG-3 pathway maybe an important immune evasion mechanism and could contribute to current cancer immunotherapy efforts for several reasons. FGL1 inhibits antigen-specific T cell activation, and ablation of FGL1 in mice promotes T cell immunity. Blockade of the FGL1-LAG-3 interaction by monoclonal antibodies stimulates tumor immunity and is therapeutic against established mouse tumors in a receptor-ligand inter-dependent manner. FGL1 is highly produced by human cancer cells, and elevated FGL1 in the plasma of cancer patients is associated with a poor prognosis and resistance to anti-PD-1/ B7-H1 therapy
基因ID	2267

基因名	FGL1
Swiss	Q08830
别名	FGL1 (YD35877)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn