

— ABLYBIO, Help Your Research



# CEACAM6 (YD35803) Rabbit mAb

货号: **AYD11651**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	37kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Mouse liver
细胞定位	Cell membrane, Apical cell membrane, Cell surface
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a protein that belongs to the carcinoembryonic antigen (CEA) family whose members are glycosyl phosphatidyl inositol (GPI) anchored cell surface glycoproteins. Members of this family play a role in cell adhesion and are widely used as tumor markers in serum immunoassay determinations of carcinoma. This gene affects the sensitivity of tumor cells to adenovirus infection. The protein encoded by this gene acts as a receptor for adherent-invasive E. coli adhesion to the surface of ileal epithelial cells in patients with Crohn's disease. This gene is clustered with genes and pseudogenes of the cell adhesion molecules subgroup of the CEA family on chromosome 19. [provided by RefSeq, Apr 2014]
基因ID	4680
基因名	CEACAM6
Swiss	P40199 ( <a href="https://www.uniprot.org/uniprotkb/P40199/entry">https://www.uniprot.org/uniprotkb/P40199/entry</a> )
别名	CEACAM6 (YD35803),CEACAM6 (YD35803) Rabbit mAb,CEACAM6,Carcinoembryonic antigen-related cell adhesion molecule 6,Non-specific crossreacting antigen,Normal cross-reacting antigen,NCA

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)