

— ABLYBIO, Help Your Research



# DDX21 (YD20096) Rabbit mAb

货号: **AYD11661**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	87kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,MCF7,HepG2,293T,Mouse spleen,Mouse thymus
细胞定位	Nucleus, nucleolus, nucleoplasm, Cytoplasm, cytosol, Mitochondrion
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is an antigen recognized by autoimmune antibodies from a patient with watermelon stomach disease. This protein unwinds double-stranded RNA, folds single-stranded RNA, and may play important roles in ribosomal RNA biogenesis, RNA editing, RNA transport, and general transcription.
基因ID	9188
基因名	DDX21
Swiss	Q9NR30 ( <a href="https://www.uniprot.org/uniprotkb/Q9NR30/entry">https://www.uniprot.org/uniprotkb/Q9NR30/entry</a> )
别名	DDX21 (YD20096),DDX21 (YD20096) Rabbit mAb,DDX21,DEAD box protein 21,Gu-alpha,Nucleolar RNA helicase Gu,Nucleolar RNA helicase II,RH II/Gu

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)