

DDX21 (YD20096) Rabbit mAb

货号: **AYD11661**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC
推荐浓度	
理论分子量	87kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,MCF7,HepG2,293T,Mouse spleen,Mouse thymus
细胞定位	Nucleus, nucleolus, nucleoplasm, Cytoplasm, cytosol, Mitochondrion
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	DEAD box proteins, characterized by the conserved motif Asp-Glu-Ala-Asp (DEAD), are putative RNA helicases. They are implicated in a number of cellular processes involving alteration of RNA secondary structure such as translation initiation, nuclear and mitochondrial splicing, and ribosome and spliceosome assembly. Based on their distribution patterns, some members of this family are believed to be involved in embryogenesis, spermatogenesis, and cellular growth and division. This gene encodes a DEAD box protein, which is an antigen recognized by autoimmune antibodies from a patient with watermelon stomach disease. This protein unwinds double-stranded RNA, folds single-stranded RNA, and may play important roles in ribosomal RNA biogenesis, RNA editing, RNA transport, and general transcription.
基因ID	9188

基因名	DDX21
Swiss	Q9NR30
别名	DDX21 (YD20096)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn