

— ABLYBIO, Help Your Research



Phospho-VAV1 (Tyr174) (YD35675) Rabbit mAb

货号: **AYD11696**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	98kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene is a member of the VAV gene family. The VAV proteins are guanine nucleotide exchange factors (GEFs) for Rho family GTPases that activate pathways leading to actin cytoskeletal rearrangements and transcriptional alterations. The encoded protein is important in hematopoiesis, playing a role in T-cell and B-cell development and activation. The encoded protein has been identified as the specific binding partner of Nef proteins from HIV-1. Coexpression and binding of these partners initiates profound morphological changes, cytoskeletal rearrangements and the JNK/SAPK signaling cascade, leading to increased levels of viral transcription and replication. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	7409
基因名	VAV1
Swiss	P15498 (https://www.uniprot.org/uniprotkb/P15498/entry)
别名	Phospho-VAV1 (Tyr174) (YD35675),Phospho-VAV1 (Tyr174) (YD35675) Rabbit mAb,VAV1,VAV

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)