

— ABLYBIO, Help Your Research



HLTF (YD20018) Rabbit mAb

货号: **AYD11755**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	114kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,293T,LO2,A-549,HeLa,HT-1080
细胞定位	Nucleus, Chromosome
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the SWI/SNF family. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein contains a RING finger DNA binding motif. Two transcript variants encoding the same protein have been found for this gene. However, use of an alternative translation start site produces an isoform that is truncated at the N-terminus compared to the full-length protein.
基因ID	6596
基因名	HLTF
Swiss	Q14527 (https://www.uniprot.org/uniprotkb/Q14527/entry)
别名	HLTF (YD20018),HLTF (YD20018) Rabbit mAb,HLTF,DNA-binding protein/plasminogen activator inhibitor 1 regulator,HIP116,Helicase-like transcription factor,RING finger protein 80,RING-type E3 ubiquitin transferase HLTF

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)