

— ABLYBIO, Help Your Research



CIN85 (YD35226) Rabbit mAb

货号: **AYD11775**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	73kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm, cytoskeleton, Cytoplasmic vesicle membrane, Synapse, synaptosome, Cell junction, focal adhesion
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Adapter protein involved in regulating diverse signal transduction pathways. Involved in the regulation of endocytosis and lysosomal degradation of ligand-induced receptor tyrosine kinases, including EGFR and MET/hepatocyte growth factor receptor, through an association with CBL and endophilins. The association with CBL, and thus the receptor internalization, may be inhibited by an interaction with PDCD6IP and/or S PRY2. Involved in regulation of ligand-dependent endocytosis of the IgE receptor. Attenuates phosphatidylinositol 3-kinase activity by interaction with its regulatory subunit (By similarity). May be involved in regulation of cell adhesion; promotes the interaction between TTK2B and PDCD6IP. May be involved in the regulation of cellular stress response via the MAPK pathways through its interaction with MAP3K4. Is involved in modulation of tumor necrosis factor mediated apoptosis. Plays a role in the regulation of cell morphology and cytoskeletal organization. Required in the control of cell shape and migration. Has an essential role in the stimulation of B cell activation (PubMed:29636373)
基因ID	30011
基因名	SH3KBP1
Swiss	Q96B97 (https://www.uniprot.org/uniprotkb/Q96B97/entry)
别名	CIN85 (YD35226),CIN85 (YD35226) Rabbit mAb,SH3KBP1,CD2-binding protein 3,Cbl-interacting protein of 85 kDa,Human Src family kinase-binding protein 1,CIN85

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)