

CD55 (YD15347) Rabbit mAb

货号: **AYD11809**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	41kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NCI-H460,Mouse testis,Rat lung
细胞定位	Cell membrane, Secreted
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a glycoprotein involved in the regulation of the complement cascade. Binding of the encoded protein to complement proteins accelerates their decay, thereby disrupting the cascade and preventing damage to host cells. Antigens present on this protein constitute the Cromer blood group system (CROM). Alternative splicing results in multiple transcript variants. The predominant transcript variant encodes a membrane-bound protein, but alternatively spliced transcripts may produce soluble proteins.
基因ID	1604
基因名	CD55
Swiss	P08174

别名	CD55 (YD15347)
----	----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn