

CD81 (YD15729) Rabbit mAb

货号: **AYD11810**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	IHC-P ICC/IF FC
推荐浓度	
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat brain
细胞定位	Cell membrane, Basolateral cell membrane
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	<p>The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the te traspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins. This protein appears to promote muscle cell fusion and support myotube maintenance. Also it may be involved in signal transduction. This gene is localized in the tumor-suppressor gene region and thus it is a candidate gene for malignancies. Two transcript variants encoding different isoforms have been found for this gene.</p>
基因ID	975

基因名	CD81
Swiss	P60033
别名	CD81 (YD15729)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn