

— ABLYBIO, Help Your Research



# MMP14 (YD11413) Rabbit mAb

货号: **AYD11833**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	66kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-1080,Mouse lung,Mouse uterus,Rat liver
细胞定位	Cell membrane, Melanosome, Cytoplasm
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the protein encoded by this gene is a member of the membrane-type MMP (MT-MMP) subfamily; each member of this subfamily contains a potential transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. This protein activates MMP2 protein, and this activity may be involved in tumor invasion.
基因ID	4323
基因名	MMP14
Swiss	P50281 ( <a href="https://www.uniprot.org/uniprotkb/P50281/entry">https://www.uniprot.org/uniprotkb/P50281/entry</a> )
别名	MMP14 (YD11413),MMP14 (YD11413) Rabbit mAb,MMP14,MMP-X1,Membrane-type matrix metalloproteinase 1,Membrane-type-1 matrix metalloproteinase

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)