

— ABLYBIO, Help Your Research



RAP1A/B (YD20337) Rabbit mAb

货号: **AYD11842**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,HeLa,Mouse kidney
细胞定位	Cell membrane, Cytoplasm, cytosol, Cell junction
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the RAS-like small GTP-binding protein superfamily. Members of this family regulate multiple cellular processes including cell adhesion and growth and differentiation. This protein localizes to cellular membranes and has been shown to regulate integrin-mediated cell signaling. This protein also plays a role in regulating outside-in signaling in platelets. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 3, 5, 6 and 9.
基因ID	5908
基因名	RAP1B
Swiss	P61224 (https://www.uniprot.org/uniprotkb/P61224/entry)
别名	RAP1A/B (YD20337),RAP1A/B (YD20337) Rabbit mAb,RAP1B,GTP-binding protein smg p21B

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)