

— ABLYBIO, Help Your Research



# Podoplanin (YD16505) Rabbit mAb

货号: **AYD11848**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung,Rat kidney
细胞定位	Membrane, Cell projection, lamellipodium membrane, filopodium membrane, microvillus membrane, ruffle membrane, Membrane raft, Apical cell membrane, Basolateral cell membrane, invadopodium, Cytoplasm, cytosol
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor. The specific function of this protein has not been determined but it has been proposed as a marker of lung injury. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	10630
基因名	PDPN
Swiss	Q86YL7 ( <a href="https://www.uniprot.org/uniprotkb/Q86YL7/entry">https://www.uniprot.org/uniprotkb/Q86YL7/entry</a> )
别名	Podoplanin (YD16505), Podoplanin (YD16505) Rabbit mAb, PDPN, Aggrus, Glycoprotein 36, PA2.26 antigen, T1-alpha, GP36

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)