

— ABLYBIO, Help Your Research



PBK (YD16175) Rabbit mAb

货号: **AYD11860**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC IP
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,THP-1,NCI-H460
细胞定位	
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a serine/threonine protein kinase related to the dual specific mitogen-activated protein kinase kinase (MAPKK) family. Evidence suggests that mitotic phosphorylation is required for its catalytic activity. The encoded protein may be involved in the activation of lymphoid cells and support testicular functions, with a suggested role in the process of spermatogenesis. Overexpression of this gene has been implicated in tumorigenesis. Alternative splicing results in multiple transcript variants.
基因ID	55872
基因名	PBK
Swiss	Q96KB5 (https://www.uniprot.org/uniprotkb/Q96KB5/entry)
别名	PBK (YD16175),PBK (YD16175) Rabbit mAb,PBK,Cancer/testis antigen 84,MAPKK-like protein kinase,Nori-3,PDZ-binding kinase,Spermatogenesis-related protein kinase,T-LAK cell-originated protein kinase,TOPK

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)