

— ABLYBIO, Help Your Research



# Activin A Receptor Type IB (YD13681) Rabbit mAb

货号: **AYD11864**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	57kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,U-87MG,Mouse brain,Rat brain
细胞定位	Cell membrane
纯化	亲和纯化




## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes an activin A type IB receptor. Activins are dimeric growth and differentiation factors which belong to the transforming growth factor-beta (TGF-beta) superfamily of structurally related signaling proteins. Activins signal through a heteromeric complex of receptor serine kinases which include at least two type I and two type II receptors. This protein is a type I receptor which is essential for signaling. Mutations in this gene are associated with pituitary tumors. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jun 2010]
基因ID	91
基因名	ACVR1B
Swiss	P36896 ( <a href="https://www.uniprot.org/uniprotkb/P36896/entry">https://www.uniprot.org/uniprotkb/P36896/entry</a> )
别名	Activin A Receptor Type IB (YD13681),Activin A Receptor Type IB (YD13681) Rabbit mAb,ACVR1B,Activin receptor type IB,Activin receptor-like kinase 4,Serine/threonine-protein kinase receptor R2,ACVRLK4,ALK4

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)