

— ABLYBIO, Help Your Research



OX40L (YD19592) Rabbit mAb

货号: **AYD11869**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,PC-3,HepG2,293T,Mouse thymus,Mouse spleen
细胞定位	Membrane
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a cytokine of the tumor necrosis factor (TNF) ligand family. The encoded protein functions in T cell antigen-presenting cell (APC) interactions and mediates adhesion of activated T cells to endothelial cells. Polymorphisms in this gene have been associated with Sjogren's syndrome and systemic lupus erythematosus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]
基因ID	7292
基因名	TNFSF4
Swiss	P23510 (https://www.uniprot.org/uniprotkb/P23510/entry)
别名	OX40L (YD19592),OX40L (YD19592) Rabbit mAb,TNFSF4,Glycoprotein Gp34,OX40 ligand,TAX transcriptionally-activated glycoprotein 1,TXGP1

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)