

— ABLYBIO, Help Your Research



Elongation Factor 1A2 (YD20077) Rabbit mAb

货号: **AYD11870**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,MCF7,293T,NCI-H460,Mouse brain,Mouse heart,Mouse skeletal muscle,Rat heart,Rat skeletal muscle
细胞定位	Endoplasmic reticulum membrane
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes an isoform of the alpha subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This isoform (alpha 2) is expressed in brain, heart and skeletal muscle, and the other isoform (alpha 1) is expressed in brain, placenta, lung, liver, kidney, and pancreas. This gene may be critical in the development of ovarian cancer.
基因ID	1917
基因名	EEF1A2
Swiss	Q05639 (https://www.uniprot.org/uniprotkb/Q05639/entry)
别名	Elongation Factor 1A2 (YD20077),Elongation Factor 1A2 (YD20077) Rabbit mAb,EEF1A2,Eukaryotic elongation factor 1 A-2,Statin-S1,EEF1AL,STN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)