

— ABLYBIO, Help Your Research



GATA1 (YD11404) Rabbit mAb

货号: **AYD11898**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	Nucleus
纯化	亲和纯化

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a protein which belongs to the GATA family of transcription factors. The protein plays an important role in erythroid development by regulating the switch of fetal hemoglobin to adult hemoglobin. Mutations in this gene have been associated with X-linked dyserythropoietic anemia and thrombocytopenia.
------	--

基因ID	2623
基因名	GATA1
Swiss	P15976 (https://www.uniprot.org/uniprotkb/P15976/entry)
别名	GATA1 (YD11404),GATA1 (YD11404) Rabbit mAb,GATA1,Eryf1,GATA-binding factor 1,NF-E1 DNA-binding protein,GF1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)