

— ABLYBIO, Help Your Research



Aurora B (YD13156) Rabbit mAb

货号: **AYD11909**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	39kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,NIH/3T3,C6
细胞定位	Nucleus, Chromosome, centromere, kinetochore, Cytoplasm, cytoskeleton, spindle, Midbody
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the aurora kinase subfamily of serine/threonine kinases. The genes encoding the other two members of this subfamily are located on chromosomes 19 and 20. These kinases participate in the regulation of alignment and segregation of chromosomes during mitosis and meiosis through association with microtubules. A pseudogene of this gene is located on chromosome 8. Alternatively spliced transcript variants have been found for this gene.
基因ID	9212
基因名	AURKB
Swiss	Q96GD4 (https://www.uniprot.org/uniprotkb/Q96GD4/entry)
别名	Aurora B (YD13156),AIK2,AIM-1,AIM1,ARK-2,ARK2,AurB,AURKB,AURKB_HUMAN,Aurora 1,Aurora and Ipl1 like midbody associated protein 1,Aurora kinase B,Aurora related kinase 2,Aurora- and Ipl1-like midbody-associated protein 1,Aurora-B,Aurora-related kinase 2,Aurora/IPL1 related kinase 2,Aurora/IPL1-related kinase 2,IPL1,PPP1R48,Protein phosphatase 1 regulatory subunit 48,Serine/threonine kinase 12,Serine/threonine protein kinase 12,Serine/threonine-protein kinase 12,Serine/threonine-protein kinase aurora-B,STK-1,STK1,STK12,STK5,Aurora B (YD13156) Rabbit mAb,Serine/threonine-protein kinase 5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)