

PKD1 (YD13185) Rabbit mAb

货号: **AYD11914**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	102kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	Cytoplasm, Cell membrane, Golgi apparatus, trans-Golgi network
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a serine/threonine protein kinase involved in many cellular processes , including Golgi body membrane integrity and transport, cell migration and differentiation, MAPK8/JNK1 and Ras pathway signaling, MAPK1/3 (ERK1/2) pathway signaling, cell survival, and regulation of cell shape and adhesion.
基因ID	5587
基因名	PRKD1
Swiss	Q15139

别名	PKD1 (YD13185)
----	----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn