

— ABLYBIO, Help Your Research



Collagen I (YD32808) Rabbit mAb

货号: **AYD11923**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	139kDa/129kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted, extracellular space, extracellular matrix
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Type I collagen is a member of group I collagen (fibrillar forming collagen) Type I collagen is a member of group I collagen (fibrillar forming collagen)
基因ID	1277, 1278

基因名	COL1A1, COL1A2
Swiss	P02452 (https://www.uniprot.org/uniprotkb/P02452/entry), P08123 (https://www.uniprot.org/uniprotkb/P08123/entry)
别名	Collagen I (YD32808), Collagen I (YD32808) Rabbit mAb, COL1A1, COL1A2, Alpha-1 type I collagen, Alpha-2 type I collagen

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)