

— ABLYBIO, Help Your Research



HDAC1 (YD13518) Rabbit mAb

货号: **AYD11944**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF FC IP
推荐浓度	
理论分子量	55kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa,HepG2,MCF7,Mouse thymus,Mouse testis,Rat thymus,Rat testis
细胞定位	Nucleus
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	Histone acetylation and deacetylation, catalyzed by multisubunit complexes, play a key role in the regulation of eukaryotic gene expression. The protein encoded by this gene belongs to the histone deacetylase/acuc/apha family and is a component of the histone deacetylase complex. It also interacts with retinoblastoma tumor-suppressor protein and this complex is a key element in the control of cell proliferation and differentiation. Together with metastasis-associated protein-2, it deacetylates p53 and modulates its effect on cell growth and apoptosis.
基因ID	3065
基因名	HDAC1
Swiss	Q13547 (https://www.uniprot.org/uniprotkb/Q13547/entry)
别名	HDAC1 (YD13518),HDAC1 (YD13518) Rabbit mAb,HDAC1,Protein deacetylase HDAC1,Protein deacylase HDAC1,RPD3L1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)