

— ABLYBIO, Help Your Research



# ABCG1 (YD13569) Rabbit mAb

货号: **AYD12044**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	76kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2
细胞定位	Endoplasmic reticulum membrane, Golgi apparatus membrane, Cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the White subfamily. It is involved in macrophage cholesterol and phospholipids transport, and may regulate cellular lipid homeostasis in other cell types. Six alternative splice variants have been identified.
基因ID	9619
基因名	ABCG1
Swiss	P45844 ( <a href="https://www.uniprot.org/uniprotkb/P45844/entry">https://www.uniprot.org/uniprotkb/P45844/entry</a> )
别名	ABCG1 (YD13569), ABCG1 (YD13569) Rabbit mAb, ABCG1, ATP-binding cassette transporter 8, White protein homolog, ABC8, WHT1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)