

— ABLYBIO, Help Your Research



# LIPH (YD35421) Rabbit mAb

货号: **AYD12080**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	51kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,DU145,293T,BxPC-3,A-549,Mouse testis,Mouse lung,Rat testis
细胞定位	Secreted, Cell membrane
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a membrane-bound member of the mammalian triglyceride lipase family. It catalyzes the production of 2-acyl lysophosphatidic acid (LPA), which is a lipid mediator with diverse biological properties that include platelet aggregation, smooth muscle contraction, and stimulation of cell proliferation and motility.
------	--

基因ID	200879
基因名	LIPH
Swiss	Q8WWY8 ( <a href="https://www.uniprot.org/uniprotkb/Q8WWY8/entry">https://www.uniprot.org/uniprotkb/Q8WWY8/entry</a> )
别名	LIPH (YD35421),LIPH (YD35421) Rabbit mAb,LIPH,LPD lipase-related protein,Membrane-associated phosphatidic acid-selective phospholipase A1-alpha,Phospholipase A1 member B,LPDLR,MPAPLA1,PLA1B

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)