

— ABLYBIO, Help Your Research



# Integrin alpha V (YD34607) Rabbit mAb

货号: **AYD12090**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	116kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,NIH/3T3,C6
细胞定位	Cell membrane, Cell junction, focal adhesion
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The product of this gene belongs to the integrin alpha chain family. Integrins are heterodimeric integral membrane proteins composed of an alpha subunit and a beta subunit that function in cell surface adhesion and signaling. The encoded preproprotein is proteolytically processed to generate light and heavy chains that comprise the alpha V subunit. This subunit associates with beta 1, beta 3, beta 5, beta 6 and beta 8 subunits. The heterodimer consisting of alpha V and beta 3 subunits is also known as the vitronectin receptor. This integrin may regulate angiogenesis and cancer progression. Alternative splicing results in multiple transcript variants. Note that the integrin alpha 5 and integrin alpha V subunits are encoded by distinct genes.
基因ID	3685
基因名	ITGAV
Swiss	P06756 ( <a href="https://www.uniprot.org/uniprotkb/P06756/entry">https://www.uniprot.org/uniprotkb/P06756/entry</a> )
别名	Integrin alpha V (YD34607), Integrin alpha V (YD34607) Rabbit mAb, ITGAV, Vitronectin receptor, Vitronectin receptor subunit alpha, MSK8

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)