

ID4 (YD18279) Rabbit mAb

货号: **AYD12097**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Nucleus
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the inhibitor of DNA binding (ID) protein family. These proteins are basic helix-loop-helix transcription factors which can act as tumor suppressors but lack DNA binding activity. Consequently, the activity of the encoded protein depends on the protein binding partner.
基因ID	3400
基因名	ID4
Swiss	P47928
别名	ID4 (YD18279)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn