

# Spermine synthase (YD19135) Rabbit mAb

货号: **AYD12119**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF
推荐浓度	
理论分子量	41kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,293T,A375,Raji,Mouse brain,Mouse heart,Mouse kidney
细胞定位	
纯化	

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a protein belonging to the spermidine/spermin synthase family. Pseudogenes of this gene are located on chromosomes 1, 5, 6 and X. Mutations in this gene are associated with X-linked Snyder-Robinson mental retardation syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	6611
基因名	SMS
Swiss	P52788

别名	Spermine synthase (YD19135)
----	-----------------------------

### 产品验证



### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)