

SATB2 (YD11970) Rabbit mAb

货号: **AYD12124**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	83kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2
细胞定位	Nucleus matrix
纯化	

抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a DNA binding protein that specifically binds nuclear matrix attachment regions. The encoded protein is involved in transcription regulation and chromatin remodeling. Defects in this gene are associated with isolated cleft palate and cognitive disability. Alternate splicing results in multiple transcript variants that encode the same protein.
基因ID	23314
基因名	SATB2
Swiss	Q9UPW6

别名	SATB2 (YD11970)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn