

— ABLYBIO, Help Your Research



Phospho-c-Jun (Ser63) (YD13611) Rabbit mAb

货号: **AYD12132**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	36kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HT-1080,NIH/3T3
细胞定位	Nucleus
纯化	亲和纯化

抗原信息




抗原信息	
------	--

靶点信息

研究背景	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
------	---

基因ID	3725
基因名	JUN
Swiss	P05412 (https://www.uniprot.org/uniprotkb/P05412/entry)
别名	Phospho-c-Jun (Ser63) (YD13611),Phospho-c-Jun (Ser63) (YD13611) Rabbit mAb,JUN,Activator protein 1,P roto-oncogene c-Jun,Transcription factor AP-1 subunit Jun,V-jun avian sarcoma virus 17 oncogene homolo g,p39

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)