

— ABLYBIO, Help Your Research



GAP43 (YD34707) Rabbit mAb

货号: **AYD12305**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Neuro-2a,Mouse heart,Rat brain
细胞定位	Cell membrane, Cell projection, growth cone membrane, Synapse, filopodium membrane, Perikaryon, dendrite, axon, Cytoplasm
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene has been termed a 'growth' or 'plasticity' protein because it is expressed at high levels in neuronal growth cones during development and axonal regeneration. This protein is considered a crucial component of an effective regenerative response in the nervous system. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
基因ID	2596
基因名	GAP43
Swiss	P17677 (https://www.uniprot.org/uniprotkb/P17677/entry)
别名	GAP43 (YD34707),GAP43 (YD34707) Rabbit mAb,GAP43,Axonal membrane protein GAP-43,Growth-associated protein 43,Neural phosphoprotein B-50,pp46

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)