

— ABLYBIO, Help Your Research



# DROSHA (YD33325) Rabbit mAb

货号: **AYD12329**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	159kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Nucleus, nucleolus, Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a ribonuclease (RNase) III double-stranded RNA-specific ribonuclease and subunit of the microprocessor protein complex, which catalyzes the initial processing step of microRNA (miRNA) synthesis. The encoded protein cleaves the stem loop structure from the primary microRNA (pri-miRNA) in the nucleus, yielding the precursor miRNA (pre-miRNA), which is then exported to the cytoplasm for further processing. In a human cell line lacking a functional copy of this gene, canonical miRNA synthesis is reduced. Somatic mutations in this gene have been observed in human patients with kidney cancer.
基因ID	29102
基因名	DROSHA
Swiss	Q9NRR4 ( <a href="https://www.uniprot.org/uniprotkb/Q9NRR4/entry">https://www.uniprot.org/uniprotkb/Q9NRR4/entry</a> )
别名	DROSHA (YD33325),DROSHA (YD33325) Rabbit mAb,DROSHA,Protein Drosha,Ribonuclease III,p241,RN3,RNASE3L,RNASEN

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)