

— ABLYBIO, Help Your Research



hnRNP Q (YD14900) Rabbit mAb

货号: **AYD12361**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF
推荐浓度	
理论分子量	70kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,A-431,MCF7,BT-474,Mouse heart,Mouse brain,Rat heart
细胞定位	Cytoplasm, Microsome, Endoplasmic reticulum, Nucleus, nucleoplasm
纯化	亲和纯化




抗原信息


抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the cellular heterogeneous nuclear ribonucleoprotein (hnRNP) family. hn RNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA) and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. The encoded protein plays a role in multiple aspects of mRNA maturation and is associated with several multiprotein complexes including the apoB RNA editing-complex and survival of motor neurons (SMN) complex. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the short arm of chromosome 20.
基因ID	10492
基因名	SYNCRIP
Swiss	O60506 (https://www.uniprot.org/uniprotkb/O60506/entry)
别名	hnRNP Q (YD14900),hnRNP Q (YD14900) Rabbit mAb,SYNCRIP,Glycine- and tyrosine-rich RNA-binding protein,NS1-associated protein 1,Synaptotagmin-binding,cytoplasmic RNA-interacting protein,HNRPQ,NSAP1

产品验证





实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)