

CD7 (YD16763) Rabbit mAb

货号: **AYD12382**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC IP
推荐浓度	
理论分子量	25kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,U-937,Rat spleen
细胞定位	Membrane
纯化	


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.
基因ID	924
基因名	CD7
Swiss	P09564
别名	CD7 (YD16763)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn