

— ABLYBIO, Help Your Research



CD62L/L-selectin (YD33290) Rabbit mAb

货号: **AYD12383**

产品信息

| | |
|-------|---|
| 反应 | Mouse |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB FC |
| 推荐浓度 | |
| 理论分子量 | 42kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | |
| 细胞定位 | Cell membrane |
| 纯化 | 亲和纯化 |

抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|------|--|
| 研究背景 | Calcium-dependent lectin that mediates cell adhesion by binding to glycoproteins on neighboring cells. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes (PubMed:1693096). Promotes initial tethering and rolling of leukocytes in endothelia (By similarity) |
|------|--|

| | |
|-------|---|
| 基因ID | 3134 |
| 基因名 | Sell |
| Swiss | P18337 (https://www.uniprot.org/uniprotkb/P18337/entry) |
| 别名 | CD62L/L-selectin (YD33290),CD62L/L-selectin (YD33290) Rabbit mAb,Sell,CD62 antigen-like family member L,Leukocyte adhesion molecule 1,Leukocyte-endothelial cell adhesion molecule 1,Lymph node homing receptor,Lymphocyte antigen 22,Lymphocyte surface MEL-14 antigen,Lnhr,Ly-22,Ly22 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)