

CD63 (YD14178) Rabbit mAb

货号: **AYD12387**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	26kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A375
细胞定位	Cell membrane, Lysosome membrane, Late endosome membrane, Endosome, multivesicular body, Melanosome, Secreted, extracellular exosome, Cell surface
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. The encoded protein is a cell surface glycoprotein that is known to complex with integrins. It may function as a blood platelet activation marker. Deficiency of this protein is associated with Hermansky-Pudlak syndrome. Also this gene has been associated with tumor progression. Alternative splicing results in multiple transcript variants encoding different protein isoforms.
基因ID	967

基因名	CD63
Swiss	P08962
别名	CD63 (YD14178)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn