

— ABLYBIO, Help Your Research



# Interferon Regulatory Factor 4 (YD31552) Rabbit mAb

货号: **AYD12395**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,Mouse lung
细胞定位	Nucleus, Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene belongs to the IRF (interferon regulatory factor) family of transcription factors, characterized by a unique tryptophan pentad repeat DNA-binding domain. The IRFs are important in the regulation of interferons in response to infection by virus, and in the regulation of interferon-inducible genes. This family member is lymphocyte specific and negatively regulates Toll-like-receptor (TLR) signaling that is central to the activation of innate and adaptive immune systems. A chromosomal translocation involving this gene and the IgH locus, t(6;14)(p25;q32), may be a cause of multiple myeloma. Alternatively spliced transcript variants have been found for this gene.
基因ID	3662
基因名	IRF4
Swiss	Q15306 ( <a href="https://www.uniprot.org/uniprotkb/Q15306/entry">https://www.uniprot.org/uniprotkb/Q15306/entry</a> )
别名	Interferon Regulatory Factor 4 (YD31552), Interferon Regulatory Factor 4 (YD31552) Rabbit mAb, IRF4, Lymphocyte-specific interferon regulatory factor, Multiple myeloma oncogene 1, NF-EM5, MUM1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)