

— ABLYBIO, Help Your Research



Calmodulin (YD16867) Rabbit mAb

货号: **AYD12405**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P FC IP
推荐浓度	
理论分子量	17kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCI-H460,U-251MG,Mouse lung,Mouse brain,Rat lung
细胞定位	Cytoplasm, cytoskeleton, spindle, spindle pole, microtubule organizing center, centrosome, Cell projection, cilium, flagellum
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a member of the EF-hand calcium-binding protein family. It is one of three genes which encode an identical calcium binding protein which is one of the four subunits of phosphorylase kinase. Two pseudogenes have been identified on chromosome 7 and X. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	801
基因名	CALM1
Swiss	P0DP23 (https://www.uniprot.org/uniprotkb/P0DP23/entry)
别名	Calmodulin (YD16867), Calmodulin (YD16867) Rabbit mAb, CALM1, CALM, CAM, CAM1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)