

Nova1 (YD14738) Rabbit mAb

货号: **AYD12408**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	52kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,SH-SY5Y
细胞定位	Nucleus
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene encodes a neuron-specific RNA-binding protein, a member of the Nova family of paraneoplastic disease antigens, that is recognized and inhibited by paraneoplastic antibodies. These antibodies are found in the sera of patients with paraneoplastic opsoclonus-ataxia, breast cancer, and small cell lung cancer. Alternatively spliced transcripts encoding distinct isoforms have been described.
基因ID	4857
基因名	NOVA1
Swiss	P51513

别名	Nova1 (YD14738)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn