

— ABLYBIO, Help Your Research



# Phospho-delta 1 Catenin (Tyr228) (YD12730) Rabbit mAb

货号: AYD12434

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	108kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,NIH/3T3,C6
细胞定位	Cell junction, adherens junction, Cytoplasm, Nucleus, Cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene.
基因ID	1500
基因名	CTNND1
Swiss	O60716 ( <a href="https://www.uniprot.org/uniprotkb/O60716/entry">https://www.uniprot.org/uniprotkb/O60716/entry</a> )
别名	Phospho-delta 1 Catenin (Tyr228) (YD12730), Phospho-delta 1 Catenin (Tyr228) (YD12730) Rabbit mAb, CTNND1, Cadherin-associated Src substrate, p120 catenin, p120(cas), KIAA0384

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)