

— ABLYBIO, Help Your Research



# Fibronectin (YD13473) Rabbit mAb

货号: **AYD12443**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	272kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,NIH/3T3
细胞定位	Secreted, extracellular space, extracellular matrix
纯化	亲和纯化

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes fibronectin, a glycoprotein present in a soluble dimeric form in plasma, and in a dimeric or multimeric form at the cell surface and in extracellular matrix. The encoded preproprotein is proteolytically processed to generate the mature protein. Fibronectin is involved in cell adhesion and migration processes including embryogenesis, wound healing, blood coagulation, host defense, and metastasis. The gene has three regions subject to alternative splicing, with the potential to produce 20 different transcript variants, at least one of which encodes an isoform that undergoes proteolytic processing. The full-length nature of some variants has not been determined.
基因ID	2335
基因名	FN1
Swiss	P02751 ( <a href="https://www.uniprot.org/uniprotkb/P02751/entry">https://www.uniprot.org/uniprotkb/P02751/entry</a> )
别名	Fibronectin (YD13473),Fibronectin (YD13473) Rabbit mAb,FN1,Cold-insoluble globulin,FN

## 产品验证

Human plasma



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)