

— ABLYBIO, Help Your Research



Cleaved N-terminal DFNA5/GSDME (YD19870) Rabbit mAb

货号: AYD12444

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	55kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	Cell membrane, Cytoplasm, cytosol
纯化	亲和纯化

抗原信息

抗原信息	
------	--

靶点信息

研究背景	Hearing impairment is a heterogeneous condition with over 40 loci described. The protein encoded by this gene is expressed in fetal cochlea, however, its function is not known. Nonsyndromic hearing impairment is associated with a mutation in this gene. Three transcript variants encoding two different isoforms have been found for this gene.
基因ID	1687
基因名	GSDME
Swiss	O60443 (https://www.uniprot.org/uniprotkb/O60443/entry)
别名	Cleaved N-terminal DFNA5/GSDME (YD19870),Cleaved N-terminal DFNA5/GSDME (YD19870) Rabbit mAb, GSDME,Inversely correlated with estrogen receptor expression 1,Non-syndromic hearing impairment protein 5,DFNA5

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn//www.ablybio.cn>)