

— ABLYBIO, Help Your Research



Dnmt3a (YD19643) Rabbit mAb

货号: **AYD12509**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	102kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,293T,HeLa,HL-60,MCF7,HepG2,Mouse testis
细胞定位	Nucleus, Chromosome, Cytoplasm
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation. Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a DNA methyltransferase that is thought to function in de novo methylation, rather than maintenance methylation. The protein localizes to the cytoplasm and nucleus and its expression is developmentally regulated.
基因ID	1788
基因名	DNMT3A
Swiss	Q9Y6K1 (https://www.uniprot.org/uniprotkb/Q9Y6K1/entry)
别名	Dnmt3a (YD19643), Dnmt3a (YD19643) Rabbit mAb, DNMT3A, Cysteine methyltransferase DNMT3A, DNA methyltransferase HsaIIIA

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)